

RAPID TESTING LOG

Agency:

Location:

Determine Tests and Controls		Other Test/Manufacturer:	Other Test/Manufacturer:
Device Lot Number:		Device Lot Number:	
Device Expiration Date:		Device Expiration Date:	
Package Insert Revision Date:		Package Insert Revision Date:	
Control Lot No.		Control Lot No.	
Control Exp. Date		Control Exp. Date	

Date	Testing I.D. sticker or +/- Control	Staff Initials	Temperature	Start Time	Read Time	Internal control valid?	Result *Pos/Neg/Inv (If Pos: Ag+, Ab+, or Ag/Ab+)	Confirmatory Result** pos/neg	Comments -Indicate reason for running control -If test is invalid, indicate next steps -If rapid test is reactive, indicate whether client received confirmatory test results

*Indicate whether positive result is Ag+, Ab+, or Ag/Ab+ **Indicate whether confirmatory positive result is Acute/Early Infection or Established Infection

Date	Testing I.D. sticker or +/- Control	Staff Initials	Temperature	Start Time	Read Time	Internal control valid?	Result *Pos/Neg/Inv (If Pos: Ag+, Ab+, or Ag/Ab+)	Confirmatory Result** pos/neg	Comments -Indicate reason for running control -If test is invalid, indicate next steps -If rapid test is reactive, indicate whether client received confirmatory test results

*Indicate whether positive result is Ag+, Ab+, or Ag/Ab+ **Indicate whether confirmatory positive result is Acute/Early Infection or Established Infection

Instructions for Rapid Testing Log

Agency	Fill in name of Agency
Location	Fill in location of testing: e.g. Smith Clinic or outreach
For each rapid test fill out...	
Other Test/Manufacturer	Fill in the type of test (RST or HCV) and the manufacturer of the test
Device Lot Number and Expiration Date	Fill in Lot Number (on outside of box for Determine) and Expiration Date of test devices. For Determine, the expiration date is determined by the shortest dated test kit component, on box or pouch.
Package Insert Revision Date	Fill in the revision date of the package insert for this box of tests. Revision date is typically listed at the end of the package insert.
Control Lot No. and Expiration Date:	Fill in Lot Number (on box for Determine) and Expiration Date of most recent control performed.
Date	Fill in date of rapid test
Test ID Sticker or +/- Control	<ul style="list-style-type: none"> • If testing a Positive Control, fill in "+ Control" • If testing a Negative Control, fill in "– Control" • If testing a client sample, fill in Test ID number (or use sticker).
Staff Initials	Fill in the initials of staff conducting the rapid test
Temperature	Fill in the current temperature of the testing site
Start Time	Write the exact time that the buffer was added to the test strip for Determine or other rapid test was started.
Read Time	Write the exact time that the result was read.
Internal Control Valid	<ul style="list-style-type: none"> • If control line is present, write "Y" for yes. • If there is no control line, write "N" for no – the test is invalid. (Explain your next steps under the comment section).
Result	<ul style="list-style-type: none"> • For a non-reactive result – write "neg " • For a reactive result—write "react" • For an Ag reactive Determine result – write "Ag+" • For an Ab reactive Determine result – write "Ab+" • For an Ag/Ab reactive Determine result – write "Ag/Ab+" • For an invalid– write "inv"
Confirmatory Sample sent?	<ul style="list-style-type: none"> • If yes – write "Y" • If no – write "N" • If not applicable (in the case of controls) – write "NA"
Confirmatory Result?	<ul style="list-style-type: none"> • For a positive final result - write "pos" and indicate whether the positive result is an early (Ag only) or established (Ab positive) infection • For a negative final result – write "neg" • If not applicable (in the case of controls) – write NA
Comments	<ul style="list-style-type: none"> • If conducting either a Positive or Negative Control, indicate reason. • If test is invalid, indicate next steps • If rapid test is reactive, indicate whether client received confirmatory test results and/or next steps.